

PROCEDURE




- 2** Rotate the lid of sample extraction buffer bottle.
CAUTION: Open it away from your face and be careful not to spill any of the liquid.




- 3** Squeeze all the extraction buffer out of the bottle into the extraction tube. **CAUTION:** Avoid touching the bottle against the tube.



- 4** Find the swab in the sealed wrapper. Identify the soft, fabric tip on the swab.



- 5** Peel open the wrapper of the swab and carefully pull out the swab. **CAUTION:** Do not touch the soft, fabric tip of the swab with your hands.



- 6** Carefully insert the swab into one nostril. The swab tip should be inserted at least 2.5 cm from the edge of the nostril. Roll the swab 3-4 times along the mucosa inside the nostril to ensure that both mucus and cells are collected. Leave swab in the nostril for several seconds. Using the same swab, repeat this process for the other nostril.

CAUTION: This may feel uncomfortable. Do not insert the swab any deeper if you feel strong resistance or pain.




- 7** Place the swab with the sample into the extraction tube. Rotate the swab 3-5 times. **Leave the swab in the extraction solution for 1 minute.**



- 8** Pinch the extraction tube together with fingers while removing the swab to leave as much solution in the tube as possible. Place the swab back into the swab wrapper.



- 9** Install the nozzle cap onto the sample extraction tube tightly. Replace the tube in the tube holder.



- 10** Open the foil pouch and remove the test card. Place the test card on a flat and level surface. **CAUTION:** Once opened, the test card must be used immediately.



- 11** Invert the extraction tube and add 3 drops only of test specimen into the specimen well (S), by gently squeezing the extraction tube.
CAUTION:
Do not add drops to the larger well.
Do not add more than 3 drops
Do not agitate the tube. Avoid the formation of air bubbles.



- 12** Read the results at 15-20 minutes. **CAUTION:** Results after 20 minutes may not be accurate.